

SPECIFICATION SHEET

Product Name: free-PSA Antigen **Catalog #:** L2C00101

Description: **Free** Prostate Specific Antigen (free-PSA)

Source: Human Seminal Fluid

Format: Purified, Liquid

Concentration: Activity concentration: 1.5mg/ml (Canag EIA)

Lot #: L180801

Buffer: PBS, pH 7.4 ± 0.2

Preservative: none

Application: Specific methodologies have not been tested using this product.

Storage: Store at -20 °C. Avoid multiple freeze/thaw cycles.

All donor units tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.