

SPECIFICATION SHEET

Product Name:	McAb to Anti-HER2/neu (p185)	Catalog #:	L1C03102				
Description:	Mouse Anti Human nerve membrane protein (HER2/neu (p185)) Monoclonal Antibody						
Specificity:	Recognizes an epitope on the HER2/neu Antigen.						
Clone:	HR002						
Host Animal:	Mouse. Hybridization of X63-Ag8 myeloma cells with spleen cells from Balb/c mice						
Isotype:	IgG1						
Source:	Mouse Ascites						
Immunogen:	HER2/neu (P185)						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration&Lot:	5.4mg/ml	Lot:	L180701				
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA,CLIA,RIA etc.Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C03102</td><td>L1C03101</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C03102	L1C03101
<u>Capture</u>	<u>Detection</u>						
L1C03102	L1C03101						
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY.