

SPECIFICATION SHEET

Product Name:	Human AFP Antigen Purified	Catalog #:	L2C00201
Description:	Human Alpha Fetoprotein (AFP) Antigen Purified AFP is a glycoprotein containing about 4% carbohydrate. AFP is produced by the fetal liver and is not normally expressed in adults.		
Source:	Human fetal cord serum		
Purification:	90% pure (SDS-PAGE). ion-exchange, affinity and gel-filtration chromatography.		
Concentration:	Activity: 1mg/ml (EIA)	Lot #:	L190501
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Immunogen, calibrator in immunoassays, in vitro studies		
Storage:	Upon receipt, store at -20°C.		
Warnings:	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HbsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		