

SPECIFICATION SHEET

Catalog #:	L2C004	Lot #:	L180701
Description:	hCG Antigen Human Chorionic Gonadotropin (hCG Antigen)		
Source:	Human Pregnancy Urine		
Format:	Purified, Liquid		
Purification:	≥80% pure (SDS-PAGE)		
Concentration:	2mg/ml (about 5,000 IU/mg)		
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4,		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		