
SPECIFICATION SHEET

Product Name:	CEAAntigen Calibrator Grade	Catalog #:	L2C010
Description:	Carcinoembryonic Antigen (CEA)		
Source:	Human colon cancer		
Format:	Liquid		
Purification:	General purification.		
Concentration:	1.0 mg/ml	Lot:	L171201
Buffer:	0.01M PBS PH7.4		
Preservative:	0.05% NaN ₃		
Applications:	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	Upon receipt, store at -20°C.		
Inactivation:	Not applicable.		
Warnings:	This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.