

**SPECIFICATION SHEET**

<b>Description:</b>	Human Liver Ferritin	<b>Catalog #:</b>	L2C013
<b>Source:</b>	Human Liver		
<b>Format:</b>	Purified, Liquid,		
<b>Concentration:</b>	2mg/ml (Roche assay)	<b>Lot:</b>	L171201
<b>Buffer:</b>	0.01M PBS PH7.4, 0.05% NaN <sub>3</sub>		
<b>Preservative:</b>	0.05% NaN <sub>3</sub>		
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage:</b>	Store at -20°C for at least three years.		
<b>Inactivation:</b>	Not applicable		
<b>Warnings:</b>	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg and HIV-1 Antigen and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**