

SPECIFICATION SHEET

Product Name:	McAb-sST2	Catalog #:	L3C01201
Description:	McAb to sST2 antibody Mouse Anti Soluble ST2 Polyclonal Antibody,		
Specificity:	Recognizes an epitope on the sST2 Antigen.		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human sST2Antigen		
Format:	Affinity Purified, Liquid,Protein A chromatography		
Purification:	≥90% pure (SDS-PAGE).		
Concentration:	4.38mg/ml (OD280nm)	Lot :	L180601
Buffer:	0.01M PB, 150mm NaCl,pH7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	store at -20°C.		