

## SPECIFICATION SHEET

<b>Product Name:</b>	PcAb to IL6	<b>Catalog #:</b>	L3C01303
<b>Description:</b>	Rabbit Anti Human IL6 Antibody,		
<b>Specificity:</b>	Recognizes an epitope on the IL6 Antigen.		
<b>Clone:</b>	None		
<b>Host Animal:</b>	Rabbit		
<b>Isotype:</b>	IgG		
<b>Source:</b>	serum		
<b>Immunogen:</b>	Recombinant IL6 Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	≥90% pure (SDS-PAGE).		
<b>Concentration&amp;Lot:</b>	7.5 mg/ml	<b>Lot :</b>	L170401
<b>Buffer:</b>	0.01M PBS, pH 7.4 ± 0.2		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application.		
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		