

SPECIFICATION SHEET

| | | | |
|--------------------------------|--|-------------------|----------|
| Product Name: | McAb to SAA | Catalog #: | L3C02004 |
| Description: | Mouse Anti Serum amyloid A (SAA) Antibody, | | |
| Specificity: | Recognizes an epitope on Serum amyloid A (SAA) Antigen. | | |
| Clone: | SA004 | | |
| Host Animal: | Mouse | | |
| Isotype: | IgG | | |
| Source: | serum | | |
| Immunogen: | Recombinant Serum amyloid A (SAA) Antigen | | |
| Format: | Purified, Liquid | | |
| Purification: | ≥90% pure (SDS-PAGE). | | |
| Concentration&Lot: | 5 mg/ml | Lot : | L190301 |
| Buffer: | 0.01M PBS, pH 7.4 ± 0.2 | | |
| Preservative: | 0.01% NaN3 | | |
| Applications: | Suitable for use in CLIA. Each laboratory should determine an optimum working titer for use in its particular application. | | |
| Storage&Shelf life: | Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles. | | |