
SPECIFICATION SHEET

Catalog #:	L4C00201	Lot #:	L180703
Description:	Troponin T-Cardiac		
Source:	Heart tissue		
Format:	Liquid		
Concentration:	0.15mg/ml		
Buffer:	20mM Tris buffer, pH 7.5 (5mM EDTA, 15mM beta-mercaptoethanol, 0.1%SKL).		
Preservative:	None		
Applications:	Suitable for use as a standard in immunoassay, immunogen for antisera production, tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage:	Store at $< -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles.		
Inactivation:	Not applicable		
Warnings:	A blood sample from the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.