

SPECIFICATION SHEET

Product Name:	Recombinant D-dimmer Antigen	Catalog #:	L4C00802
Description:	Recombinant Human D-Dimer protein		
Source:	Expressed in E. coli		
Purification:	Purified, >90% pure by SDS-PAGE		
Concentration&Lot:	1.0 mg/ml	Lot:	L181001
Buffer:	50mM Tris-Hcl, 300mM NaCl, 0.1%SKL ,PH 8.0		
Preservative:	0.1%SKL		
Applications:	Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.