

SPECIFICATION SHEET

Product Name:	PCT Antigen	Catalog #:	L4C00901
Description:	PCT Antigen, Recombinant procalcitonin, (PCT) protein (16KD)		
Source:	E.coli		
Purification:	Affinity purified.		
Concentration&Lot:	0.65mg/ml (Roche ECLIA)	Lot:	L190201
Buffer:	20MM Citric Acid, pH 4.5		
Preservative:	None		
Applications:	Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.