

SPECIFICATION SHEET

**Catalog #:** L4I01101

**Description:** Recombinant RBP-4 ( Retinol-Binding Protein 4,RBP 4)

**SYNONYMS:** Retinol Binding Protein 4, RBP-4, RBP4, Plasma retinol-binding protein, PRBP, RBP.

**Source:** E. Coli **Species:** Homo sapiens

**Molecular weight:** 21kd (with His Tag protein)

**Concentration:** 1mg/ml

**Lot #:** L170801

**Purity :** >95% by SDS-PAGE

**Formulation :** The Protein solution contains 20mM saline sodium citrate pH3.5, 1mm EDTA,5%Trehalose

**Preservative:** None

**Application:** ELISA; Positive Control; Immunogen; SDS-PAGE; WB. Other applications have not been tested . If bio-activity of the protein is needed, please check active protein.

**Storage & Expiration:** Upon receipt, aliquot and store at <-20°C .for long term.Avoid repeated freeze and thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.