

SPECIFICATION SHEET

Product Name:	McAb to HE4	Catalog #:	L1C02601				
Description:	Mouse Anti Human epididymis protein 4 (HE4) Monoclonal Antibody						
Specificity:	Recognizes an epitope on the HE4 Antigen.						
Clone:	HE001						
Host Animal:	Mouse. Hybridization of X63-Ag8 myeloma cells with spleen cells from Balb/c mice						
Isotype:	IgG						
Source:	Mouse Ascites						
Immunogen:	Human HE4 Antigen						
Format:	Purified, Liquid						
Purification:	Protein A chromatography						
Concentration&Lot:	6.5mg/ml	Lot :	L180401				
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA, CLIA, RIA etc. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C02601</td><td>L1C02602</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C02601	L1C02602
<u>Capture</u>	<u>Detection</u>						
L1C02601	L1C02602						
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY.