

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb to Pro-GRP	<b>Catalog #:</b>	L1C03202				
<b>Description:</b>	Mouse Anti Human Pro-GRP Monoclonal Antibody						
<b>Specificity:</b>	Recognizes an epitope on the Human Pro-GRP Antigen.						
<b>Clone:</b>	GR002						
<b>Host Animal:</b>	Mouse. Hybridization of X63-Ag8 myeloma cells with spleen cells from Balb/c mice						
<b>Isotype:</b>	IgG						
<b>Source:</b>	Mouse Ascites						
<b>Immunogen:</b>	Human Pro-GRP Antigen						
<b>Format:</b>	Purified, Liquid						
<b>Purification:</b>	Protein A chromatography						
<b>Concentration&amp;Lot:</b>	4.2 mg/ml	<b>Lot:</b>	L180201				
<b>Buffer:</b>	0.01M PBS, 0.15M NaCl, pH 7.4						
<b>Preservative:</b>	None						
<b>Applications:</b>	Suitable for use in ELISA, CLIA, RIA etc. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C03201</td><td>L1C03202</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C03201	L1C03202
<u>Capture</u>	<u>Detection</u>						
L1C03201	L1C03202						
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.						

**FOR RESEARCH USE ONLY.**