

SPECIFICATION SHEET

Product Name:	McAb to MCM5	Catalog #:	L1C03301				
Description:	McAb to Human MCM5 (Human Minichromosome Maintenance Deficient 5) Monoclonal Antibody						
Specificity:	Human Minichromosome Maintenance Deficient 5. No cross-reactivity with Human Albumin, Other cytokines.						
Clone:	MC001						
Host Animal:	Mouse. Hybridization of P3X63-Ag8.653 .myeloma cells with spleen cells from Balb/c mice.						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Purified MCM5 from Cord Blood						
Format:	Purified, Liquid						
Purification:	DEAE-chromatography						
Concentration&Lot:	3.5mg/ml	Lot #:	L190701				
Buffer:	0.01M PB, 0.15M NaCl, pH 7.4.						
Preservative:	None						
Applications:	Determination of MCM5. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C03301</td><td>L1C03302</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C03301	L1C03302
<u>Capture</u>	<u>Detection</u>						
L1C03301	L1C03302						
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.						