

SPECIFICATION SHEET

Description: Recombinant human HE4 Antigen **Catalog #:** L2C02601

Molecular Weight: 14kd (calculated from sequence)
14KD(determined by SDS gelelectrophoresis)

Source: E. coli

Format: Purified, Liquid

Purification: > 90% (determined by SDS gelelectrophoresis)

Concentration: 1mg/ml (OD280 OD260) **Lot:** L160701

Buffer: 20Mm Tris-HCL, pH 8.0 + 0.2

Preservative: 0.02%NaN3

Storage: Store at -20 C for 3years. Avoid multiple freeze/thaw cycles.