

SPECIFICATION SHEET

Product Name:	HBE-Ag	Catalog #:	L1H00501
Description:	Recombinant Hepatitis B e Ag Hepatitis B e Antigen (Hepatitis B “e” Antigen, 160aa,17kD)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE)		
Concentration:	1mg/ml	Lot.:	L180501
Buffer:	0.02M Sodium phosphate, pH 7.2		
Preservative:	0.02% Sodium azide		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years to 2020-1. Avoid multiple freeze/thaw cycles.		
Warnings:	All human source materials have tested negative for HIV1, HIV2 and HCV antibodies. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

For Research Use Only.