

SPECIFICATION SHEET

<b>Product Name:</b>	McAb to HBsAg	<b>Catalog #:</b>	L1H00202				
<b>Description:</b>	Monoclonal Antibody to Hepatitis B Surface Antigen (HBsAg)						
<b>Specificity:</b>	HBsAg.						
<b>Clone:</b>	N/A						
<b>Host Animal:</b>	Mouse. Hybridization of Sp2/0 myeloma						
<b>Isotype:</b>	IgG cells with spleen cells from BALB/c mice.						
<b>Source:</b>	Ascites						
<b>Immunogen:</b>	Native HBsAg						
<b>Format:</b>	Purified, Liquid						
<b>Purification:</b>	>90% pure. Protein A chromatography						
<b>Concentration:</b>	2.4mg/ml	<b>Lot:</b>	L190702				
<b>Buffer:</b>	PBS, pH 7.4						
<b>Preservative:</b>	None						
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1H00202</td><td>L1H00201</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1H00202	L1H00201
<u>Capture</u>	<u>Detection</u>						
L1H00202	L1H00201						

**Storage& Shelf life:** Store at -20° C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.