

SPECIFICATION SHEET

Product Name: McAb to HBcAg **Catalog #:** L1H00601

Description: Monoclonal Antibody to Hepatitis B c Antigen (HBcAg)

Specificity: HBcAg.

Clone: LT05

Host Animal: Mouse. Hybridization of Sp2/0 myeloma

Isotype: IgG cells with spleen cells from BALB/c mice.

Source: Ascites

Immunogen: HBcAg

Format: Purified, Liquid

Purification: >90% pure. Protein A chromatography

Concentration: 2.5mg/ml **Lot.:** H141001

Buffer: 10MM PBS, pH 7.4

Preservative: None

Applications: Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage& Shelf life: Store at -20° C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.