

SPECIFICATION SHEET

Product Name:	TRecombinant TB-16&38	Catalog #:	L1H01302
Description:	Recombinant TB-16 and TB-38 antigen, Antigen Grade		
Source:	E. coli		
Format:	Purified, Protein solution in Tris-HCl		
Purification:	Purity >90%		
Concentration&Lot:	5mg/ml	Lot:	L161201
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.02,		
Preservative:	None		
Applications:	EIA,CLIA,IP,ICH,WB		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		