

## SPECIFICATION SHEET

<b>Product Name:</b>	Recombinant HP UreaseA	<b>Catalog #:</b>	L1H01203
<b>Description:</b>	Rcombinant Helicobacter pylori (Hp) Urease A Antigen, Antigen Grade		
<b>Molecular Weight:</b>	26.71KD (calculated from sequence) 26.7KD (determined by SDS gelelectrophoresis)		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	> 85% (determined by SDS gelelectrophoresis)		
<b>Concentration&amp;Lot:</b>	1.29mg/ml,	<b>Lot:</b>	L180401
<b>Buffer:</b>	20Mm Tris-HCl, pH 8.0 ± 0.1		
<b>Preservative:</b>	0.2%NaN3		
<b>Application:</b>	Immunoassays or Immunization.		
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		