

SPECIFICATION SHEET

Product Name:	Syphilis TP17 antigen	Catalog #:	L2H00701
Description:	Recombinant syphilis TP17 antigen,Antigen Grade		
Source:	E. coli		
Format:	Purified, Protein solution in Tris-HCl pH8.0		
Purification:	Purity >90%, 17KD		
Concentration&Lot:	1~5mg/ml, lot specific		
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.02,		
Preservative:	None		
Application:	Immunoassays as EIA, IFA or Immunization.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		