

SPECIFICATION SHEET

Product Name:	TP62(47&15)	Catalog #:	L2H00703
Description:	Recombinant syphilis TP-47and TP-15 antigen, Antigen Grade		
Source:	E. coli		
Format:	Purified, Protein solution in Tris-Glycin pH8.0		
Purification:	Purity >90%		
Concentration&Lot:	1~5mg/ml, lot specific		
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.02,		
Preservative:	None		
Application:	Immunoassays as EIA, IFA or Immunization.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		