

SPECIFICATION SHEET

Product Name: EV71 VP1 Recombinant Protein **Catalog #:** L9H001

SYNONYMS: P1D, Virion protein 1

Source: E. Coli **Species:** Human enterovirus 71

Molecular weight: 37KD (with His Tag protein)

Concentration: 2.55mg/ml **Lot #:** L181201

Purity : >85% by SDS-PAGE

Formulation : The Protein solution contains 20mMTrisHCl,PH=8.0, 1mm EDTA,5%Trehalose,25mM L-Arg

Application: ELISA; Positive Control; Immunogen; SDS-PAGE; WB. Other applications have not been tested .

Storage & Expiration: Upon receipt, aliquot and store at <-20°C .for long term.Avoid repeated freeze and thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.