

SPECIFICATION SHEET

Product Name:	McAb to NT-Pro-BNP	Catalog #:	L3C00504
Description:	Monoclonal Antibody to Human NT-pro-Brain Natriuretic Peptide (NT-proBNP)		
Specificity:	NT-proBNP		
Clone:	NB004		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Recombinant NT-pro-Brain Natriuretic Peptide (NT-proBNP)		
Format:	Purified, Liquid		
Purification:	Protein G chromatography		
Concentration:	5 mg/ml	Lot:	L190101
Buffer:	10mM PBS, pH 7.4		
Preservative:	none		
Applications:	Can be used in an immunoassay for specific determination of human myoglobin. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay:		
Storage:	Short term store at 2-8°C. Long term store at -20°C. Avoid multiple/freeze thaw cycles.		