

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb to NT-Pro-BNP	<b>Catalog #:</b>	L3C00506
<b>Description:</b>	Monoclonal Antibody to Human NT-pro-Brain Natriuretic Peptide (NT-proBNP)		
<b>Specificity:</b>	NT-proBNP		
<b>Clone:</b>	NB006		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Recombinant NT-pro-Brain Natriuretic Peptide (NT-proBNP)		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	Protein G chromatography		
<b>Concentration&amp;Lot:</b>	7.5mg/ml	<b>Lot :</b>	L181001
<b>Buffer:</b>	10mM PBS, pH 7.4		
<b>Preservative:</b>	none		
<b>Applications:</b>	Can be used in an immunoassay for specific determination of human myoglobin. Each laboratory should determine an optimum working titer for use in its particular application.		
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		