

SPECIFICATION SHEET

Product Name:	McAb-FABP	Catalog #:	L3C00801
Description:	McAb to FABP antibody Mouse Anti Human Fatty Acid Binding Protein (FABP) Monoclonal Antibody,		
Specificity:	Recognizes an epitope on the FABP Antigen.		
Clone:	FA001		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human FABP Antigen		
Format:	Affinity Purified, Liquid, Protein A chromatography		
Purification:	≥90% pure (SDS-PAGE).		
Concentration:	2.0mg/ml	Lot:	L181201
Buffer:	0.01M PB, 150mm NaCl, pH7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY.