

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb-FABP	<b>Catalog #:</b>	L3C00802
<b>Description:</b>	McAb to FABP antibody Mouse Anti Human Fatty Acid Binding Protein (FABP) Monoclonal Antibody,		
<b>Specificity:</b>	Recognizes an epitope on the FABP Antigen.		
<b>Clone:</b>	FA002		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Purified Human FABP Antigen		
<b>Format:</b>	Affinity Purified, Liquid, Protein A chromatography		
<b>Purification:</b>	≥90% pure (SDS-PAGE).		
<b>Concentration:</b>	2.7mg/ml	<b>Lot:</b>	L170501
<b>Buffer:</b>	0.01M PB, 150mm NaCl, pH7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		

**Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.