

SPECIFICATION SHEET

Product Name:	McAb-LpPLA2	Catalog #:	L3C01001
Description:	McAb to LpPLA2 antibody Mouse Anti Human Phospholipase A2, Lipoprotein Associated (LpPLA2) Monoclonal Antibody,		
Specificity:	Recognizes an epitope on the LpPLA2 Antigen.		
Clone:	LP001		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human LpPLA2 Antigen		
Format:	Affinity Purified, Liquid, Protein A chromatography		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	5 mg/ml	Lot:	L180501
Buffer:	0.01M PB, 150mm NaCl, pH7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY.