

SPECIFICATION SHEET

Product Name:	McAb-PCT	Catalog #:	L3C00904
Description:	McAb to PCT antibody Mouse Anti Human procalcitonin (PCT) Monoclonal Antibody,		
Specificity:	Recognizes an epitope on the PCT Antigen.		
Clone:	PC004		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human PCT Antigen		
Format:	Affinity Purified, Liquid, Protein A chromatography		
Purification:	≥90% pure (SDS-PAGE).		
Concentration:	1.25mg/ml	Lot:	L160101
Buffer:	0.01M PB, 150mm NaCl, pH7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.