

SPECIFICATION SHEET

| | | | |
|--------------------------------|--|-------------------|----------|
| Product Name: | McAb to IL6 | Catalog #: | L3C01301 |
| Description: | Mouse Anti Human IL6 Antibody, | | |
| Specificity: | Recognizes an epitope on the IL6 Antigen. | | |
| Clone: | I6001 | | |
| Host Animal: | Mouse | | |
| Isotype: | IgG | | |
| Source: | serum | | |
| Immunogen: | Recombinant IL6 Antigen | | |
| Format: | Purified, Liquid | | |
| Purification: | ≥90% pure (SDS-PAGE). | | |
| Concentration&Lot: | 5.2 mg/ml | Lot : | L181001 |
| Buffer: | 0.01M CPBS, pH 6.0 ± 0.2 | | |
| Preservative: | 0.01% NaN3 | | |
| Applications: | Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. | | |
| Storage&Shelf life: | Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles. | | |