

SPECIFICATION SHEET

Product Name:	McAb to IL6	Catalog #:	L3C01305
Description:	Mouse Anti Human IL6 Antibody,		
Specificity:	Recognizes an epitope on the IL6 Antigen.		
Clone:	I6005		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	serum		
Immunogen:	Recombinant IL6 Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	4.6mg/ml	Lot :	L190601
Buffer:	0.01M PBS, pH 6.0 ± 0.2		
Preservative:	0.01% NaN3		
Applications:	Suitable for use in LFA. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		