

SPECIFICATION SHEET

Product Name:	PcAb to SAA	Catalog #:	L3C02003
Description:	Rabbit Anti Human SAA Antibody,		
Specificity:	Recognizes an epitope on the SAA Antigen.		
Clone:	None		
Host Animal:	Rabbit		
Isotype:	IgG		
Source:	serum		
Immunogen:	Recombinant Human SAA Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	4.5 mg/ml	Lot :	L170901
Buffer:	0.01M PBS, pH 7.4 ± 0.2		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		