

Recombinant human SAA protein.

- Catalog #:** L4C01201 **Lot #:** L190501
- Description:** Serum amyloid A (SAA), the precursor protein in inflammation-associated reactive amyloidosis (AA-type)
- Source:** E. Coli
- Molecular weight:** 16KD (with Tag protein)
- Concentration:** 2mg/ml
- Purity:** >85% by SDS-PAGE
- Formulation:** Liquid
- Application:** ELISA WB Antigen Standards
- Storage buffer:** 20mM PH7.4, 150mM NaCl, 0.05% NaN₃
- Storage & Expiration:** Ship at 4°C. Upon receipt, aliquot and store at -20°C or -40°C for long term. Avoid repeated freeze and thaw cycles.

For research use only.