

SPECIFICATION SHEET

- Product Name:** Tox-Ag P30 **Catalog #:** L1Y00301
- Description:** Recombinant T. gondii P30 Ag (P30,SAG1)
- Source:** E.Coli, 30KD
Strain: RH from ATCC
- Format:** Lysate, Liquid
- Purification:** >85% (SDS-PAGE)
- Concentration&Lot.:** 2mg/ml (Lot. L161001)
- Buffer:** PBS, 0.15M NaCl, pH 7.4 ± 0.2
- Preservative:** None
- Applications:** Specific methodologies have not been tested using this product.
- Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.