

SPECIFICATION SHEET

Product Name: GAD230 **Catalog #:** L4I00303

Description: Recombinant human Glutamic Acid Decarboxylase(GAD230)

Molecular Weight: 45kd (calculated from sequence)
45KD(determined by SDS gelectrophoresis)

Source: E. coli

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelectrophoresis)

Concentration& Lot.: 1.47mg/ml (OD280) (Lot. specific)

Buffer: Tris-HCL, pH 7.4 ± 0.4

Preservative: 0.02%NaN3

Application: Immunoassays Immunization.

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY