

SPECIFICATION SHEET

Product Name: ICA 69 **Catalog #:** L4I00403

Description: Recombinant human Islet cells Antigen(ICA 69)

Source: Pichia pastoris , MW: about 72KD

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelelectrophoresis)

Concentration& Lot.: 1.66mg/ml **Lot.:** L190401

Buffer: 10mM TrisHCl,0.1%SKL,2mmEDTA

Preservative: 0.1%SKL

Application: Immunoassays Immunization.

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.