

SPECIFICATION SHEET

Description: Recombinant human TSH protein. **Catalog #:** L2G00102

Description: Recombinant Thyroid Stimulating Hormone,

Source: E. Coli

Molecular weight: 48 KD (with Tag protein)

Concentration: 1.5mg/ml **Lot #:** L180901

Purity : >85% by SDS-PAGE

Formulation : Liquid

Application: ELISA; Other applications have not been tested. The optimal dilutions should be determined by end user.

Storage buffer: 10mM PBS (pH7.4 ± 0.2)

Storage & Expiration: Upon receipt, aliquot and store at <-20°C
for long term. Avoid repeated freeze and thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.