

SPECIFICATION SHEET

Catalog #:	L2G00201	Lot #:	L170501
Description:	Follicle Stimulating Hormone (FSH) intact human		
Purity:	>85% by SDS-PAGE		
Source:	Human pituitary glands		
Concentration:	1mg/ml (1800U/mg)		
Contaminants:	not detectable		
Buffer:	10mM PBS, 150mM NaCl, PH 7.4		
Biohazard:	Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.		
Storage:	-20°C or lower, for long term storage.		