

SPECIFICATION SHEET

Product Name: McAb to Recombinant ABCA7 Protein **Catalog #:** L1B00801

Description: McAb to Recombinant ABCA7 Protein
Monoclonal Antibody to ATP Binding Cassette Subfamily A
Member 7

Specificity: Reacts with free ABCA7 Protein..

Clone: AB001

Host Animal: Mouse. Hybridization of Sp2/0 myeloma

Isotype: IgG₁, cells with spleen cells from Balb/c mice.

Source: Ascites

Immunogen: Recombinant ABCA7 Protein

Format: Purified, Liquid

Purification: Protein A chromatography

Concentration&Lot: 3.9 mg/ml **Lot:** L80401

Buffer: 10mM PBS, pH 7.4

Preservative: none

Applications: Suitable for use in immunoassay

Storage&Shelf life: Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.