

SPECIFICATION SHEET

Product Name:	Rabbit anti Recombinant ABCA7 Protein	Catalog #: L1B00803
Description:	PcAb to Recombinant ABCA7 Protein Rabbit anti-ATP Binding Cassette Subfamily A Member 7 Polyclonal antibody	
Specificity:	Reacts with Recombinant ABCA7 Protein.	
Clone:	N/A	
Host Animal:	Rabbit	
Source:	Serum	
Immunogen:	Recombinant ABCA7 Protein.	
Format:	Purified, Liquid	
Purification:	Protein A chromatography	
Concentration&Lot:	3.2 mg/ml	Lot: L190301
Buffer:	10mM PBS, pH 7.4	
Preservative:	none	
Applications:	Suitable for use in immunoassay,	
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.