

SPECIFICATION SHEET

Product name: Human AD7c-NTP **Catalog #:** L2B007

Protein description: Human AD7c-NTP (AD7thcloneneuronalthreadprotein) epitope peptide1

Molecular weight 50KD **Expression host** E.coli

Purity: > 90 % by SDS-PAGE

Biological activity A identified biomarker of Alzheimer's disease

Concentration& Lot.: 1mg/ml **Lot.:** L180601

Storage buffer 0.01M PBS, 0.15M NaCl, pH 7.4, 2.5% glycerol, 0.2%ykl

Preservative 0.2%ykl

Application ELISA ;WESTERNUsed in diagnosisand of Alzheimer's disease

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.