

SPECIFICATION SHEET

Catalog #:	L1K00101	Lot #:	L160301
Description:	McAb to Cystatin C antibody Mouse Anti Human Cystatin C Monoclonal Antibody		
Specificity:	Recognizes an epitope on Cystatin C Antigen		
Clone:	CY001		
Host Animal:	Mouse		
Source:	Ascites		
Immunogen:	Recombinant Cystatin C Antigen		
Format:	Purified, Liquid		
Purification:	Protein A chromatography		
Concentration:	3.8mg/ml		
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4		
Preservative:	None		
Applications:	Suitable for use in WB, IHC		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		