

SPECIFICATION SHEET

| | | | |
|--------------------------------|--|---------------|---------|
| Catalog #: | L1K00104 | Lot #: | L180301 |
| Description: | McAb to Cystatin C antibody Mouse Anti Human Cystatin C Monoclonal Antibody | | |
| Specificity: | Recognizes an epitope on Cystatin C Antigen | | |
| Clone: | CY004 | | |
| Host Animal: | Mouse | | |
| Source: | Ascites | | |
| Immunogen: | Recombinant Cystatin C Antigen | | |
| Format: | Purified, Liquid | | |
| Purification: | Protein A chromatography | | |
| Concentration: | 3.25mg/ml | | |
| Buffer: | 0.01M PBS, 0.15M NaCl, pH 7.4 | | |
| Preservative: | None | | |
| Storage&Shelf life: | Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles. | | |