

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb to Glycated hemoglobin	<b>Catalog #:</b>	L1K00504
<b>Description:</b>	Mouse anti-Glycated hemoglobin (hemoglobin A1c, HbA1c or HGBA1C) antibody		
<b>Specificity:</b>	Human		
<b>Clone:</b>	HB001		
<b>Host Animal:</b>	Mouse.		
<b>Isotype:</b>	IgG <sub>1</sub>		
<b>Immunogen:</b>	Purified Glycated hemoglobin (hemoglobin A1c, HbA1c or HGBA1C)		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4. No preservatives added.		
<b>Concentration&amp; Lot.:</b>	4.6mg/ml	<b>Lot. :</b>	L190301
<b>Buffer:</b>	PBS, pH 7.4		
<b>Applications:</b>	Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		