

SPECIFICATION SHEET

Product Name:	Rabbit anti Sheep IgG (H&L)	Catalog #:	L1I01001
Description:	Rabbit Antibody to Sheep IgG, Heavy & Light chains		
Specificity:	Reacts with heavy chains on mouse IgG and with the light chains on all mouse immunoglobulins.		
Host Animal:	Rabbit		
Immunogen:	Purified Sheep IgG		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE). Purified.		
Concentration:	9mg/ml	Lot.:	L180301
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
Preservative:	None		
Applications:	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY