

## SPECIFICATION SHEET

**Product Name:** Rabbit IgG **Catalog #:** L11018

**Description:** Rabbit ImmunoglobulinG(IgG)

**Isotype:** IgG

**Source:** Normal Rabbit Serum..

**Immunogen:** Purified Rabbit IgG

**Format:** 0.01M PBS pH 7.4

**Concentration:** 10mg/ml (OD280nm) **Lot.:** L170601

**Purification:** Octoic acid purified.

**Application:** Sepecific methodologies have not been tested using this product.

**Preservative:** None

**Storage& Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

**FOR RESEARCH USE ONLY**