

SPECIFICATION SHEET

Product Name: Rabbit anti Folic acid (Folate)-BSA **Catalog #:** L6B01103

Synonyms: Rabbit anti-BSA conjugated FA

Specificity: Recognizes an epitope on the Folic acid (Folate)-BSA

Immunogen: Folic acid (Folate)-BSA

Host Animal: Rabbit

Source: Serum

Format: Affinity Purified, Liquid

Concentration: 5.4mg/ml **Lot:** L181101

Buffer: 0.01M PBS, pH 7.4

Preservative: None

Description: Folate, distinct forms of which are known as folic acid, folacin, and vitamin B₉, is one of the B vitamins. Synthetic vitamin B₉. Tetrahydrofolate precursor. Promotes DNA synthesis and cell division. Shows anti-inflammatory and chemopreventive effects *in vivo*.

Storage&Shelf life: Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.