

SPECIFICATION SHEET

**Catalog #:** L1K00203 **Lot #:** L20200301  
**Description:** Rabbit anti BSA-S-adenosyl-L-homocysteine (BSA-SAH)  
Polyclonal Antibody  
**Specificity:** Recognizes an epitope on the BSA-S-adenosyl-L-homocysteine  
**Immunogen:** BSA-S-adenosyl-L-homocysteine (BSA-SAH)  
**Host Animal:** Rabbit  
**Source:** Serum  
**Format:** Purified, Liquid  
**Concentration:** 4 mg/ml  
**Buffer:** 0.01M PBS, pH 7.4  
**Preservative:** None  
**Applications:** ELISA  
**Storage:** Store at -20 deg C for long term storage. Avoid repeated freeze-thaw cycles.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**